
HemoCue® Hb 201⁺

Operating Manual
Bedienungsanleitung
Manuel d'utilisation
Gebruiksaanwijzing

Table of Contents	Inhaltsverzeichnis	Table des matières	Inhoudsopgave	
HemoCue® Hb 201+ System	HemoCue® Hb 201+ System	HemoCue® Hb 201+ System	HemoCue® Hb 201+ System	1
Analyzer Overview	Überblick über den Analyzer	Présentation de l'analyseur	Overzicht analyser	3
Components	Bestandteile	Composants	Componenten	4
Start-up	Inbetriebnahme	Démarrage	Opstarten	5
Set-up	Einstellung	Configuration	Set-up	6
Printer and Connectivity	Drucker und Konnektivität	Imprimante et connectivité	Printer en connectiviteit	7
Specimen Collection Capillary	Probenentnahme Kapillär	Prélèvement des échantillons Capillaires	Bloedafname Capillair	8
Specimen Collection Venous/Arterial	Probenentnahme Venös/Arteriell	Prélèvement des échantillons Veineux/Artériels	Bloedafname Veneus/arterieel	10
Measurement Procedure	Messverfahren	Procédure de mesure	Meetprocedure	11
Quality Control	Qualitätskontrolle	Contrôle qualité	Kwaliteitscontrole	12
View Results	Ergebnisse anzeigen	Afficher les résultats	Resultaten bekijken	13
Delete Results	Ergebnisse löschen	Supprimer des résultats	Resultaten verwijderen	14
Maintenance Daily	Pflege Täglich	Maintenance Quotidienne	Onderhoud Dagelijks	15
Maintenance Optical parts	Pflege Optische Teile	Maintenance Pièces optiques	Onderhoud Optisch gedeelte	16
Troubleshooting	Fehlersuche	Dépannage	Problemen oplossen	17
Performance, EMC and Electrical Safety	Leistung, EMV und elektrische Sicherheit	Performances, CEM et sécurité électrique	Prestaties, EMC en elektrische veiligheid	25
Symbols used	Verwendete Symbole	Symboles utilisés	Gebruikte symbolen	29
Technical Specifications	Technische Daten	Caractéristiques techniques	Technische specificaties	30

HemoCue® Hb 201⁺ System

EN

Intended Purpose

The HemoCue® Hb 201⁺ System is intended for quantitative determination of hemoglobin in capillary, venous or arterial whole blood using a specially designed analyzer, the HemoCue Hb 201⁺ Analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201⁺ Analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

Intended User and Indications for Use

The HemoCue Hb 201⁺ System is an automated system for professional use, intended for near-patient (point-of-care) and laboratory testing.

Indications for use is quantitative determination of hemoglobin to support clinical decisions in the detection of anemias, erythrocytosis and polycythemias.

IVD Medical Device Regulation

The HemoCue Hb 201⁺ System complies with the IVD Medical Device Regulation (EU) 2017/746 and carries the CE mark.

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Warranty

The analyzer carries a 24-month warranty from the day of receipt. Any other use of the system than recommended by the manufacturer, including opening the cover of the analyzer, will void the warranty.

Service and Disposal

The analyzer should be cleaned prior to service. Consult local environmental authorities for proper disposal of the analyzer, batteries and other related components. Decontamination is required before disposal (see section *Maintenance*).

 Read and follow this operating manual and the HemoCue Hb 201 Microcuvettes package insert to attain optimum performance and safety. Any other use of the system than recommended by the manufacturer may impair the safety.

Further information may be obtained from the HemoCue distributor or HemoCue AB.

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

Materials provided

1. HemoCue Hb 201+ Analyzer
2. Power adapter
3. HemoCue Hb 201+ Operating Manual & Quick Reference Guide

Materials required but not provided

HemoCue Hb 201 Microcuvettes
Lancet*
Pipette or other transfer device**
Lint-free wipe
Hydrophobic surface**

Spare parts

Battery Lid
Power adapter***
Cuvette Holder

Optional items

Printer
Software applications
Cables for PC and printer connection
HemoCue Cleaner
Hard carrying case
HemoCue Lancet
HemoCue Hb 201+ Control Cuvette

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
Always follow operating and storage conditions listed in section *Technical Specifications*. Allow the analyzer to reach operating temperature before use.

Place the analyzer on a horizontal and stable surface. Connect the power adapter to the power inlet or insert 4 batteries, see section *Analyzer Overview*, fig. 2.

⚠ Only use power adapter and batteries listed in section *Technical Specifications*.


1. In Loading position (see also section *Analyzer Overview*, fig. 1C), press left button until all symbols appear on the display.

2. The software version is displayed followed by time & date, audio signal "█", an hourglass "█", and "Hb". During this time the analyzer performs a self-test.

3. Three flashing dashes indicate that the analyzer has passed the self-test and is ready for use.

To turn the analyzer off press left button until OFF is displayed. When turned on again the most recent result is displayed if the cuvette holder is in Measuring position (see section *Analyzer Overview*, fig. 3).

Low battery is indicated with a battery symbol "█".

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

Set-up / Einstellung / Configuration / Set-up

- In Loading position (see section *Analyzer Overview*, fig. 1C), enter set-up by pressing both buttons until a flashing QC symbol is displayed.
- Audio signal (e.g. to indicate measurement is completed): press right button to scroll until the audio symbol is displayed. Use left button to a) activate or b) deactivate.
- Time & Date: press right button to scroll until characters for time, year and date are displayed. The hour character is flashing, press left button to change (hold for quick advancement). Press right button to continue to next character.
- Exit set-up by pressing right button until three flashing dashes are displayed.

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Only use ASCII printers recommended by HemoCue. Always follow instructions for use for the printer.

1. Connect analyzer to printer or computer* via serial or USB port (analyzer can be turned on or off).
2. When measurement is completed, the result will automatically be printed or transferred to the computer with date and time.

Time and date must be activated before data can be printed.

* The analyzer can be connected to HemoCue software or external software applications for transfer of data, using either a USB cable or serial cable. For more information, contact the local distributor or HemoCue AB.

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2

3

Be aware that peripheral circulatory failure of the patient could affect the result of a capillary sample. To minimize variation in capillary sampling results, carefully follow the step-by-step procedure for how to perform a capillary sample, as described in this section.

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1. Make sure that the patient's hand is warm and relaxed. Use the middle or ring finger for sampling. Avoid fingers with rings on.
2. Clean fingertip with disinfectant and allow to dry.
3. Using your thumb, lightly press the finger from the top of the knuckle towards fingertip to stimulate blood flow.

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


6

7

4. Sample at the side of the fingertip for best blood flow and comfort.
5. Press lightly towards fingertip and puncture using a lancet.
6. Wipe away the first 2 or 3 drops of blood. Press lightly towards fingertip until another drop of blood appears.
7. When the blood drop is large enough, fill the microcuvette completely in one step. **Do not refill.**

Follow section *Measurement Procedure* to perform measurement.

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

⚠ Always wear protective gloves. Handle blood with care, as it may be infectious. Follow local safety procedures for disposal of used microcuvettes.

Use anticoagulants (e.g. EDTA or heparin), preferably in solid form to avoid dilutional effects. Hemoglobin remains stable for at least 6 days in 2–8 °C (36–46 °F), provided the blood does not become contaminated.

Important: Fill the microcuvette within 3 minutes after removal from the package. Always keep the vial properly closed.

1. If refrigerated, allow sample to reach operating temperature before mixing. Mix thoroughly on a mixer for at least 2 minutes or gently invert manually 8–10 times.
2. Place a drop of blood onto a hydrophobic surface using a transfer device.
3. Fill the microcuvette completely in one step. **Do not refill.**

Follow section *Measurement Procedure* to perform measurement.

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3

4

Start measurement no later than 10 minutes after filling the microcuvette.

- Wipe off specimen from the outside of the microcuvette, make sure no specimen is drawn out from the open end. Visually inspect the microcuvette.

Important: If the microcuvette is not completely filled with blood or if there are air bubbles, discard and fill a new microcuvette. Small bubbles around the edges can be ignored.

- Place the microcuvette in the cuvette holder and start measurement by gently sliding the cuvette holder to Measuring position (see section *Analyzer Overview*, fig. 3).

- Whilst measuring, will be displayed.

- The result is displayed as long as the cuvette holder is in Measuring position (see section *Analyzer Overview*, fig. 3).

Discard the microcuvette after measurement. Microcuvettes are for single use only.

If not in use, an analyzer operated on batteries will automatically turn off after 5 minutes.

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

The analyzer has an internal quality control, the self-test. It automatically verifies the performance of the analyzer every time it is turned on and every second hour when in use.

If in addition an external quality control is required by local or other regulation, only use controls recommended by HemoCue.

Perform External Quality Control

Always follow instructions for use for control solution.

1. In Loading position (see section *Analyzer Overview*, fig. 1C): press both buttons until a flashing QC symbol is displayed, then left button to select QC-test (press right button to unselect).
2. Place a drop of control solution on a hydrophobic surface. Fill the microcuvette completely in one step. **Do not refill.**

Follow section *Measurement Procedure* to perform measurement. After measurement the analyzer automatically exits QC-test.

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Up to 600 results can be stored, after which the oldest will be deleted.

In Loading or Measuring position (see section *Analyzer Overview*, fig. 1C or fig. 3) press left or right button to scroll between stored results.

1. To view sequence number, press any button until displayed. When released, the result is shown. The most recent result is always sequence number 1.
2. Exit View Results: in Loading position (see section *Analyzer Overview*, fig. 1C) wait until three flashing dashes are displayed.

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


Deleted results cannot be recovered.

1. In Loading or Measuring position (see section *Analyzer Overview*, fig. 1C or fig. 3) press both buttons until a flashing QC symbol is displayed. Press right button, a flashing waste bin is displayed.
2. Delete current result (must be done immediately after measurement): press left button until a waste bin, a warning triangle and "1" are displayed followed by "0" that confirms deletion.
3. Delete all results; press both buttons until waste bin, warning triangle and total number of stored results are displayed, followed by "000" that confirms the deletion.

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The cuvette holder should be cleaned after each day of use.

1. In Loading position (see section *Analyzer Overview*, fig. 1C) press the catch and carefully rotate the cuvette holder as far as possible to the left.
2. Remove the cuvette holder completely by pulling it off the steel pin.
3. Clean the cuvette holder with alcohol (20-70%) or mild detergent.
4. The cuvette holder must be completely dry before reattached. Allow minimum 15 minutes to dry.

The cover may be cleaned with alcohol (20-70%) or mild detergent.

Analyzer disposal

Dispose of the battery lid and cuvette holder as potentially infectious waste. Remove the batteries. Clean the analyzer cover and optical parts according to section *Maintenance*. Disinfect the analyzer external surfaces and optical parts according to local procedures. Wait for three weeks, then dispose of the analyzer according to local procedures for electronic waste.

For further information, contact HemoCue distributor or HemoCue AB.

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
1b

2

Clean optical parts when required in section *Troubleshooting*.

Follow instructions thoroughly to make sure the cleaner reaches both upper and lower cover glasses.

1. Turn off the analyzer. Follow steps 1-2 in section *Maintenance Daily* to remove the cuvette holder. Push the HemoCue Cleaner as much as possible into the opening of the analyzer. a) Move from side to side 5-10 times. b) Placed to the left, push the cleaner back and forth 5-10 times, cleaning the cover glasses. If the cleaner is stained, repeat a) and b) with a new cleaner.
2. The optical parts must be completely dry before cuvette holder is reattached. Allow minimum 15 minutes to dry.

For alternative to the HemoCue Cleaner contact the local distributor or HemoCue AB.

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Troubleshooting

If the problem remains after recommended actions below are taken, contact the HemoCue distributor or HemoCue AB. The analyzer should be cleaned prior to service, see section *Maintenance* for instructions. If any spare parts are lost or damaged, see section *Components*.

Do not open the cover of the analyzer. Warranty is void if analyzer has been opened.

* Follow relevant section in this operating manual.

Error Code	Explanation	Action
Error code is displayed	May be a temporary fault or faulty electronics /optical parts.	Turn analyzer off. Turn on after 30 seconds. If the problem continues, see specific error code below or the analyzer needs service.
E00	1. Faulty microcuvette.	1a. Check expiry date. 1b. Fill a new microcuvette and perform measurement. *
E01-E02	1. Dirty optical parts. 2. Analyzer too hot/cold. 3. Magnet missing in cuvette holder.	1. Clean optical parts. * 2. Turn analyzer off, allow to reach operating temperature before use. 3. Order new cuvette holder.
E03	1. Analyzer exposed to direct light.	1. Avoid direct light exposure.
E05, E06	1. Analyzer too hot/cold or exposed to direct light.	1a. Turn analyzer off, allow to reach operating temperature before use. 1b. Avoid direct light exposure.
E07	1. Battery power too low.	1a. Batteries need to be replaced. * 1b. Use power adapter.
E08	1a. The absorbance is too high. 1b. Faulty microcuvette or sample. 1c. Dirty or blocked optical parts.	1a and 1b. Ensure analyzer and microcuvettes are used according to operating manual and relevant package insert. 1a and 1c. Clean optical parts. *
E09-E30	1. Internal error in analyzer.	1. Clean optical parts. *
E35	1. Current measurement incorrectly stored.	1. Delete all stored results.

Troubleshooting

Measuring Deviations	Explanation	Action
HHH	1. Result exceeds measuring range.	1a. Fill a new microcuvette and perform measurement.* 1b. Confirm result with laboratory method.
No characters displayed	1. Analyzer not receiving power. 2. Faulty display.	1a. Ensure power adapter is properly connected and cable not damaged. 1b. Batteries need to be replaced.* 2. Analyzer needs service.
Erroneous characters displayed	1. Faulty display/electronics.	1. Analyzer needs service.
is displayed	1. Low battery. 2. Faulty power adapter.	1. Batteries need to be replaced.* 2. Ensure power adapter is properly connected and cable not damaged.
and "Hb" continuously displayed	1. Faulty cuvette holder/analyzer.	1. The analyzer needs service.
Unexpected patient or control results	<i>Patient sample and control solution</i> 1. Microcuvettes expired, damaged or improperly stored. 2. Microcuvette contaminated. 3. Dirty optical parts. 4. Air bubbles in filled microcuvette. 5. Faulty optical parts. 6. Incorrect sampling technique. 7. Measurement started later than 10 minutes after filling the microcuvette. <i>Control solution</i> 8. Control solution not compatible. 9. Improper handling of control solution. 10. Control solution expired or improperly stored.	<i>Patient sample and control solution</i> 1. Check expiry date and storage conditions. 2. Fill a new microcuvette and perform measurement.* 3. Clean optical parts.* 4. Fill a new microcuvette and perform measurement.* 5. The analyzer needs service. 6. Fill a new microcuvette and perform measurement.* 7. Start measurement no later than 10 minutes after filling the microcuvette. <i>Control solution</i> 8. Only use control solutions recommended by HemoCue. 9. Always follow instructions for use for control solution. 10. Check expiry date and storage conditions.

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Performance, EMC and Electrical Safety

Principle of the method and procedure

The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. Blood is drawn into the microcuvette by capillary action and an azidemethemoglobin reaction occurs. Analyzer measures transmittance at two wavelengths to compensate for turbidity. The hemoglobin level is automatically calculated and presented. The system is factory calibrated against hemiglobincyanide (HiCN) method, the international reference method for hemoglobin determination recommended by International Council for Standardization in Haematology (ICSH). No further calibration is needed.

Measuring range

The system has a measuring range of 0.5 – 25.6 g/dL (5 – 256 g/L, 0.3 – 15.9 mmol/L) and a detection limit of 0.1 g/dL (1 g/L, 0.1 mmol/L).

The system is linear within displayed range 0 – 25.6 g/dL (0 – 256 g/L, 0 – 15.9 mmol/L).

For information regarding subjects below, refer to HemoCue Hb 201 Microcuvettes Package Insert

Limitations

Reference Values

Reporting of Serious Incidents

Storage and Handling (HemoCue Hb 201 Microcuvettes)

Composition

Specific Performance Characteristics

References

- HemoCue Hb 201 Microcuvettes Package Insert
- Makarem, A. In clinical Chemistry: Principles and Techniques, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelmann, J.W., Harper and Row, Hagerstown, M.D., 1974, pp. 1125–1147
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15

EMC and Electrical Safety

The analyzer has been tested and complies with following standards:

- IEC 61010-1
- IEC/EN 61010-2-101
- IEC/EN 61326-2-6 (including applicable parts of IEC/EN 61326-1)

The analyzer has been tested for indoor use.

It is recommended to evaluate the electromagnetic environment prior to use of the device.

Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g. unshielded intentional RF sources), as these can interfere with the proper operation.

Note: It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.

Note: It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.

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Literatur

- HemoCue Hb 201 Microcuvettes Packungsbeilage
- Makarem, A. In clinical Chemistry: Principles and Techniques, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelmann, J.W., Harper and Row, Hagerstown, M.D., 1974, S. 1125–1147
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15

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