

Type	<b>C</b> SERVOTEST	<b>5+N</b> SERVOTEST	<b>9</b> SERVOTEST	<b>10<sup>SG+</sup></b> SERVOTEST	<b>10</b> SERVOTEST	<b>11</b> SERVOTEST
pH-Wert	•	•	•	•	•	•
Glucose (mg/dl)	•	•	•	•	•	•
Ascorbinsäure (mg/dl)		•				•
Keton (mg/dl)			•	•	•	•
Nitrit		•	•	•	•	•
Protein (mg/dl)	•	•	•	•	•	•
Bilirubin			•	•	•	•
Urobilinogen (mg/dl)		•	•	•	•	•
Okkultes Blut (Ery/ $\mu$ l)		•	•	•	•	•
Leukozyten (WBC/ $\mu$ l)		•		•	•	•
Spezifisches Gewicht			•	•	•	•