



ZPACK13961000 PLASTIC BOTTLE LIQ 1 L DIN 45 UN X



SPECIFICATIONS

UN-certification :	3H1/X1.4/250/
Material :	HDPE
Neck Type :	DIN 45
Height :	215,00 mm
Length :	97,00 mm
Wide :	97,00 mm

AVAILABLE ACCESSORIES

Adaptator for dispenser S40 to GL 45 : 217-9271

Key for opening : 84544.001

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343